
Exploring Ce³⁺:YAG glass ceramics by light sheet microscopy

A. Herrmann, Y. P. Hartmann, J. Burkard, M. Hofmann,
U. Brokmann, E. Rädlein

Institute of Materials Science and Engineering,
Technische Universität Ilmenau, Ilmenau, Germany


60th Ilmenau Scientific Colloquium: Living Glass Surfaces XII,
07. September 2023

- explore the crystall/glass interface by light-sheet fluorescence microscopy
incl. fluorescence spectroscopy & fluorescence lifetime measurements
- sample: Ce:YAG glass ceramics as a model for material inhomogeneity
⇒ use the fluorescence of crystall and glass phases as structural probe
- work in progress:
determine boundaries of the method / push it to its limits
⇒ probable future application in material characterization
 e.g. phase separated glasses
 photo-active glasses

Light-sheet microscopy

Light-sheet fluorescence microscope (LSFM):


- excitation: laser @ 532 nm, 25 mW (405 nm, 50 mW)
- laser beam formed into a „laser sheet“ by a cylindrical lens
⇒ only a thin slide of the sample is irradiated
- fluorescence light from the sheet is filtered and detected by SCMOS-sensor
⇒ fluorescence picture or by spectrometer
⇒ emission spectrum



Light-sheet microscopy

Light-sheet fluorescence microscope (LSFM):



- excitation: laser @ 532 nm (405 nm)
- laser beam formed into a „laser sheet“ by a cylindrical lens
⇒ only a thin slice of the sample is irradiated
- fluorescence light from the sheet is filtered and detected by SCMOS-sensor
⇒ fluorescence picture or by spectrometer
⇒ emission spectrum



Light-sheet microscopy

Light-sheet fluorescence microscope (LSFM):


- fluorescence light from the sheet is filtered and detected by SCMOS-sensor (or spectrometer)
- moving the sample: different light sheets/ stacks of light sheets
⇒ 3D picture of the sample



Ce:YAG glass ceramics

Yttrium Aluminum Garnet (YAG)

- $Y_3Al_5O_{12}$, cubic, colorless
- refractive index: 1.833
- melting temperature: 1950 °C
- chemically very stable
- very good solubility for rare earth ions
⇒ Ce^{3+} doping (yellow fluorescence)
- precipitation from a glass (previous project)
⇒ Ce:YAG glass ceramics [3]




Ce:YAG glass ceramics

Ce:YAG glass ceramics

- taken from different project
- glass: yttrium aluminosilicate
doping: 1 mol% CeO₂
- YAG crystallization at only 1040 °C [3]
- some surface crystallization of YAG after casting
 - ⇒ cubic crystals („pyramids“ at surface)
 - crystallite size 100-200 μm
 - ⇒ dendritic growth in sample volume
 - ⇒ crystals are Ce³⁺ doped (fluorescence)



surface
crystallization




Ce³⁺ fluorescence

⁵⁸Ce: electron configuration [Xe] 6s² 4f²

- 5d orbital: five broad overlapping energy levels
- 4f orbital: two narrow energy levels (‘shielding’)
- 5d-4f transitions very broad and **sensitive to host material**
- excitation/emission in most glasses: UV-blue

Ce³⁺ doped glass under UV excitation





Ce³⁺ fluorescence

⁵⁸Ce: electron configuration [Xe] 6s² 4f²

- 5d orbital: five broad overlapping energy levels
- 4f orbital: two narrow energy levels (‘shielding’)
- 5d-4f transitions very broad and **sensitive to host material**

⇒ absorption and emission shift with host material [4]






Ce³⁺ fluorescence

⁵⁸Ce: electron configuration [Xe] 6s² 4f²

- 5d orbital: five broad overlapping energy levels
- 4f orbital: two narrow energy levels (‘shielding’)
- 5d-4f transitions very broad and **sensitive to host material**

⇒ absorption and emission shift with host material [4]





Ce³⁺ fluorescence

⁵⁸Ce: electron configuration [Xe] 6s² 4f²

- 5d orbital: five broad overlapping energy levels
- 4f orbital: two narrow energy levels (‘shielding’)
- 5d-4f transitions very broad and **sensitive to host material**



⇒ absorption and emission shift with host material [4]

⇒ **changes in the emission spectra at the crystal/glass interface?**




LSFM microscopic/fluorescence pictures

- preliminary results [5]:
samples with volume crystallization
- ⇒ dendritic Ce:YAG crystals
size $\sim 100\text{-}200\ \mu\text{m}$
- excitation 532 nm:
fluorescence well detectable
⇒ 3D pictures




LSFM microscopic pictures

- new samples:
mostly surface crystallization
+ some volume crystallization
- ⇒ dendritic Ce:YAG crystals
size ~100-200 μm
- ⇒ single dendrits visible




sample surface




LSFM: spatial resolution
(luminescence picture, Ex.: 532 nm)

- new samples:
mostly surface crystallization
+ some volume crystallization
- ⇒ dendritic Ce:YAG crystals
size ~100-200 μm
- ⇒ single dendrits visible




LSFM: spatial resolution
(luminescence picture, Ex.: 532 nm)

- new samples:
mostly surface crystallization
+ some volume crystallization
- ⇒ dendritic Ce:YAG crystals
size ~100-200 μm
- ⇒ single dendrits visible
- ⇒ intensity profile




LSFM: spatial resolution
(luminescence picture, Ex.: 532 nm)

- new samples:
mostly surface crystallization
+ some volume crystallization
- ⇒ dendritic Ce:YAG crystals
size $\sim 100\text{-}200\ \mu\text{m}$
- ⇒ single dendrites visible
- ⇒ intensity profile
- ⇒ distance between dendrites $\sim 1.5\text{-}2.0\ \mu\text{m}$




LSFM emission spectra

- 532 nm laser:
low fluorescence intensity
550 nm filter cuts peak maximum
⇒ not suitable to detect peak shift




LSFM emission spectra

- 532 nm laser:
low fluorescence intensity
550 nm filter cuts peak maximum
⇒ not suitable to detect peak shift
- 405 nm laser:
much higher fluorescence intensity
good S/N ratio



LSFM emission spectra

- 532 nm laser:
low fluorescence intensity
550 nm filter cuts peak maximum
⇒ not suitable to detect peak shift
- 405 nm laser:
much higher fluorescence intensity
good S/N ratio
- measurements at different crystal positions
⇒ but: peak position constant
⇒ no peak shift at crystal/glass interface detectable



Summary & Outlook

- LSFM can be used to create 3D pictures of the precipitated crystals
- spatial resolution $\sim 1\text{-}2\ \mu\text{m}$
- fluorescence spectra could be measured with 532 and 405 nm excitation (532 nm not suitable, latter has much better S/N ratio)
- up to now no shift of the emission spectra detected at the interface crystal/glass

Outlook:

- improvement of the spatial resolution (better fitting optics for 405 nm)
- scanning the samples at higher accuracy
- fluorescence lifetime measurements at different spots in the glass ceramics
- UV laser (fluorescence of glass and YAG phases detectable)

- [1] P. Pitrone, J. Schindelin, L. Stuyvenberg et al.: OpenSPIM: an open-access light-sheet microscopy platform. Nature Methods 10 (2013) p. 598-599.
- [2] Z. Yanchun, X. Huimin, F. Zhihai: Theoretical Investigation on Mechanical and Thermal Properties of a Promising Thermal Barrier Material: $\text{Yb}_3\text{Al}_5\text{O}_{12}$, Journal of Materials Science & Technology 30 (2014) p. 631.
- [3] X. Xia: Optimization of a glass composition for synthesis of Ce^{3+} :YAG crystallites: MA thesis, TU Ilmenau, 2021.
- [4] A. Herrmann, H. A. Othman, A. A. Assadi, M. Tiegel, S. Kuhn, C. Rüssel: Spectroscopic properties of cerium-doped aluminosilicate glasses, Optical Materials Express 5 (2015) p. 720.
- [5] M. Hofmann, A. Herrmann, U. Brokmann: Light-sheet fluorescence microscopic probing of silicate materials, tm - Technisches Messen 89 (2022) p. 447-454.

Thank you for your attention!
Vielen Dank für Ihre Aufmerksamkeit!

